
# AProof®

CORONAVIRUS ANTIKÖRPERTEST

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

DETAILED INSTRUCTIONS FOR USE

INSTRUCTIONS FOR USE




1. Wash your hands with soap and water. **Tip: Warm water stimulates blood circulation.** Remove the disinfectant wipe/alcohol wipe from the packaging and disinfect the fingertip with the alcohol wipe.

2. Remove the filter paper card from the test kit and place it on a clean surface. Twist the safety cap to unlock the lancet, place the lancet on the side of the chosen fingertip and press it firmly. **Note:** Only touch the filter card at the edge.

CONTENT

 <b>1× alcohol wipe</b> CE Vivomed GmbH	 <b>2× plaster</b> CE 0123 <b>STERILE EO</b> NOBAMED Paul Danz AG <b>1× cellulose pad</b> CE Lohmann & Rauscher GmbH & Co. KG
 <b>1× filter paper card</b> PerkinElmer Health Sciences, Inc.	 <b>2× automatic lancet</b> CE 0197 <b>STERILE R</b> Vivomed GmbH
 <b>1× code card</b>	 <b>1× paper bag,</b> <b>1× prepaid return envelope</b>

REFERENCE

- IVD** In-vitro-Diagnostic
-  Instructions for use
-  For single use only
-  Store between 2 °C and 30 °C
-  Contains 1 test
-  Manufacturer
-  Expiration date
-  Recyclable Material
-  Biohazard
- CE** Conformité Européenne
- LOT** Lot number
- PZN** Pharmazentralnummer
- EAN** European Article Number
- REF** Reference number
- COMP** Components
- STERILE R** Sterilized by radiation
- STERILE EO** Sterilized by ethylene oxide

SPECIFICATIONS

The **AProof®** Coronavirus Antibody Test is an ELISA method for the detection of IgG antibodies.

During validation, 100% sensitivity was demonstrated 14 days after confirmed infection. In an examination of over 1,500 control samples, a very high specificity of 99.4% was proven.

Sensitivity describes the ability of a test to detect positive samples as positive.

Specificity describes the ability of a test to detect negative samples as negative.

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**REF** 300001-001

5. Write your address well readable on the return envelope and drop it into the nearest postbox **within 48 hours**. Shipping within Germany is free of charge. **By sending in the sample material you waive the statutory right of revocation.**

6. 

Registration + Results:  
[aproof-corona.de](http://aproof-corona.de)

Create a customer account at [www.aproof-corona.de](http://www.aproof-corona.de) (choose your own user name/password), activate your customer account via an activation link (E-mail), register the test using the **test ID from your code card**. Results will be available 1-2 working days after your sample has arrived.

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**Adversis Pharma GmbH**