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Do not reuse
Store between 2-30°C
Caution
Use by
Keep away from sunlight
Manufacturer
Version No.: 1.0
Effective Date: July 20, 2020

IVD
For in vitro diagnostic use only
Consult instructions for use
Lot number
Contains sufficient for +10 tests
Keep dry
Do not use if package is damaged
109148201

ZHENRUI COVID-19 ANTIGEN TEST CASSETTE

ENGLISH

For professional use only.
For in vitro diagnostic use only.

【 Package Insert 】

A rapid test for the qualitative detection of COVID-19 antigen in nasopharyngeal swab specimens. For professional medical institutions use only. Not for self testing

【 INTENDED USE 】

The COVID-19 Antigen Rapid Test Cassette is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

【 PRINCIPLE 】

The COVID-19 Antigen Test Cassette is a qualitative membrane strip based immunoassay for the detection of COVID-19 antigen in nasopharyngeal swab specimen. In this test procedure, anti-COVID-19-N antibody is immobilized in the test line region of the device. After a nasopharyngeal swab specimen is placed in the specimen well, it reacts with anti-COVID-19-N antibody coated particles that have been applied to the specimen pad. This mixture migrates chromatographically along the length of the test strip and interacts with the immobilized anti-COVID-19 antibody. If the specimen contains COVID-19 antigen, a colored line will appear in the test line region indicating a positive result. If the specimen does not contain COVID-19 antigen, a colored line will not appear in this region indicating a negative result. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

【 PRECAUTIONS 】

- I. For professional in vitro diagnostic use only. Do not use after expiration date.
- II. Do not eat, drink or smoke in the area where the specimens and kits are handled.
- III. Handle all specimens as if they contain infectious agents.
- IV. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- V. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- VI. Follow standard biosafety guidelines for handling and disposal of potential infective material.
- VII. Humidity and temperature can adversely affect results.

【 STORAGE AND STABILITY 】

Store as packaged in the sealed pouch at room temperature or refrigerated (4-30°C). The test is stable to the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

【 SPECIMEN COLLECTION AND PREPARATION 】

- I. **Swab** Specimen collection Only the swab provided in the kit is to be used for nasopharyngeal swab collection. To collect a nasopharyngeal swab sample, carefully insert the swab into the nostril exhibiting the most visible drainage, or the nostril that is most congested if drainage is not visible. Using gentle rotation, push the swab until resistance is met at the level of the turbinates (less than one inch into the nostril). Rotate the swab 5 times or more against the nasal wall then slowly remove from the nostril. Using the same swab, repeat sample collection in the other nostril.
- II. **COVID-19 Antigen Rapid Test (Lateral Flow)** can be applied to nasopharyngeal swab.
- III. Do not return the nasopharyngeal swab to the original paper packaging.
- IV. For best performance, direct nasopharyngeal swabs should be tested as soon as possible after collection. If immediate testing is not possible, and to maintain best performance and avoid possible contamination, it is highly recommended the nasopharyngeal swab is placed in a clean, unused plastic tube labeled with patient information, preserving sample integrity, and capped tightly at room temperature (15-30°C) for up to 1 hour prior to testing. Ensure the swab fits securely within the tube and the cap is tightly closed. If greater than 1 hour delay occurs, dispose of sample. A new sample must be collected for testing.
- V. If specimens are to be transported, they should be packed in compliance with local regulations covering the transportation of etiological agents.

【 MATERIALS 】

- I. **Materials provided:** Test devices, Extraction buffer, Extraction Tubes, Package insert Swab
- II. **Materials required but not provided:** Timer

【 DIRECTIONS FOR USE 】


- Allow the test, specimen, buffer and/or controls to reach room temperature 15-30°C (59-86°F) prior to testing.
1. Insert the swab into the extraction tube which contains 0.5 mL of the extraction reagent. Rotate the swab inside the tube using a circular motion to roll the side of the extraction tube so that the solution is expressed and reabsorbed from the swab. Leave the swab in the Reagent Tube for one minute.
 2. Pinch the extraction tube with fingers and remove the solution from the swab as far as possible. The extracted solution will be used as test specimen.
 3. Keep the extraction tube tightly closed.

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- Allow the tester and samples to equilibrate to temperature (15-30 °C or 59-86 °C) before testing.
1. Remove the test cassette from the sealed pouch.
 2. Invert the sample extraction tube and hold the sample extraction tube upright, add 3 drops (approx. 100 µl) to the sample well (S) of the test cassette, and then start the timer. See the picture below.
 3. Wait for colored lines to appear. Interpret the test results after 15 minutes. Do not read the results after 20 minutes.

【 INTERPRETATION OF RESULTS 】

Positive: Two lines appear. One line should always appear in the control line region(C), and another one apparent colored line should appear in the test line region. *NOTE: The intensity of the color in the test line regions may vary depending on the concentration of COVID-19 antibodies present in the specimen. Therefore, any shade of color in the test line region should be considered positive.

Negative: One colored line appears in the control region(C). No apparent colored line appear in the test line region.

Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【 QUALITY CONTROL 】

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

【 LIMITATIONS 】

1. This test detects both viable (live) and non-viable, SARS-CoV and COVID-19. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
2. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.
3. The performance of COVID-19 Antigen Test Cassette was evaluated using the procedures provided in this product insert only. Modifications to these procedures may alter the performance of the test.
4. False negative results may occur if a specimen is improperly collected, transported, or handled.
5. False results may occur if specimens are tested past 1 hour of collection. Specimens should be tested as quickly as possible after specimen collection.
6. Positive test results do not rule out co-infections with other pathogens.
7. Positive test results do not differentiate between SARS-CoV and COVID-19.
8. Negative test results are not intended to rule in other non-SARS viral or bacterial infections.
9. Negative results, from patients with symptom onset beyond seven days, should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management, may be performed.
10. If the differentiation of specific COVID-19 viruses and strains is needed, additional testing, in consultation with state or local public health departments, is required.

【 PERFORMANCE CHARACTERISTICS 】

Sensitivity and Specificity The COVID-19 Antigen Test Cassette (nasopharyngeal Swab Specimen) was compared with a commercial PCR (from Oriental gene)

Method	Gold Standard Reagent(PCR)		Total Result
	Result	Positive	
SARS-CoV-Antigen Rapid Test COVID-19	Positive	96	99
	Negative	4	101
Total Result		100	200

Relative Sensitivity: 96.00% (95% CI: 90.20% ~ 98.40%)
Relative Specificity: 97.00% (95% CI: 91.50% ~ 99.00%)
Accuracy: 96.5% (95% CI: 93.00% ~ 98.30%)

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Concentration	Positive/Total	Conformity
50 TCID ₅₀ / ml	180/180	100%

Cross-reactivity

The COVID-19 Antigen Test Cassette has been tested for other virus (Table below). The results showed no cross-reactivity.

No cross-reactivity	Human Rhinovirus 14	Arcanobacterium	Staphylococcus aureus subsp. aureus
Human coronavirus OC43	Human Rhinovirus 16	Candida albicans	Staphylococcus epidermidis
Coronavirus NL63	Measles	Corynebacterium	Streptococcus pneumoniae
Influenza A H1N1	Mumps	Escherichia coli	Streptococcus pyogenes
Influenza A H3N2	Parainfluenza virus 2	Moraxella catarrhalis	Streptococcus salivarius
Influenza B	Parainfluenza virus 3	Neisseria lactamica	Streptococcus sp group F
Human Rhinovirus 2	Respiratory syncytial virus	Neisseria subflava	Pseudomonas aeruginosa

Cross reactivity with following Virus or Bacteria culture with certain concentration has been studied. The results were found negative when tested with the COVID-19 Antigen Rapid Test :

Virus/Bacteria	Concentration	Results
Influenza A (H1N1)	1×10 ⁶ PFU/ml	-
Influenza A (H3N2)	1×10 ⁶ PFU/ml	-
Influenza B (Yamagata)	1×10 ⁶ PFU/ml	-
Influenza B (Victoria)	1×10 ⁶ PFU/ml	-
Adenovirus	1×10 ⁶ PFU/ml	-
Human metapneumovirus	1×10 ⁶ PFU/ml	-
Parainfluenza virus	1×10 ⁶ PFU/ml	-
Respiratory syncytial virus	1×10 ⁶ PFU/ml	-
Streptococcus pyogenes	1×10 ⁶ PFU/ml	-
Candida albicans	1×10 ⁶ PFU/ml	-
Mycoplasma pneumoniae	1×10 ⁶ PFU/ml	-
Chlamydia pneumoniae	1×10 ⁶ PFU/ml	-
Legionella pneumophila	1×10 ⁶ PFU/ml	-
Human coronavirus 229E	1×10 ⁶ PFU/ml	-
Human coronavirus OC43	1×10 ⁶ PFU/ml	-
Human coronavirus NL63	1×10 ⁶ PFU/ml	-
Human coronavirus HKU1	1×10 ⁶ PFU/ml	-
Mers (Middle East Respiratory Syndrome Coronavirus)	1×10 ⁶ PFU/ml	-

【 PRECISION 】

Intra-Assay

Within-run precision has been determined by using 15 replicates of two specimens: a negative, and an N protein (1ng/ml) as positive. The specimens were correctly identified >99% of the time.

Inter-Assay

Between-run precision has been determined by 15 independent assays on the same two specimens: a negative, an N protein (1mg/ml) as positive. Three different lots of the SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) have been tested using these specimens. The specimens were correctly identified >99% of the time.

CE IVD i 2

Shenzhen Zhenrui Biotech Co., Ltd
Add: 5F, Building B, Fuanna Industrial Park, No.1 Qinning Road, Longhua District, Shenzhen, China

CITOTEST LABWARE MANUFACTURING Co., Ltd.
No. 48, XINXU ROAD HAIMEN CITY JIANGSU PROVINCE CHINA

EC REP CE

Obelis s.a.
Bd. Général Wallez 53
B-1030 Brussels, Belgium
Phone: 32.2.732.59.54
Fax: 32.2.732.60.03
E-mail: mail@obelis.net

EC REP CE 0197

Wellkang Ltd
Enterprise Hub, NW Business Complex, 1 Beraghmore Rd Derry, BT48 8SE Northern Ireland

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Do not reuse
Store between 2-30°C
Caution
Use by
Keep away from sunlight
Manufacturer

For in vitro diagnostic use only
Consult instructions for use
LOT
Contains sufficient for $n-1$ tests
Keep dry
Do not use if package is damaged

Version No.: 1.0
Effective Date: July 20, 2020
109148201

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Importer: ALIC Deutschland GmbH